




**YOU CAN
HOME MEDICAL**



TABLE OF CONTENTS



02.	About Us
03.	Rental Equipment
04.	Oxygen & Respiratory
05.	CPAP & Sleep Therapy
06.	Walkers & Wheelchairs
07.	Lift Chairs & Recliners
09.	Mobility Scooters & Chairs
11.	Seat Cushions & Pillows
	Wellness & Pain
12.	Aids to Daily Living
13.	Compression Stockings
14.	Orthopedic & Recovery
15.	Home Safety
17.	Urology & Incontinence
	Catheters
18.	Hospital Beds

About You Can Home Medical



You Can Home Medical








in Denver, Colorado, we treat our customers like family. From walkers, wheelchairs, oxygen, compression and more. You Can Home Medical has been locally owned and operated since 1995. We are the highest rated home medical supply store in the Denver Metro area. We pride ourselves on the absolute best in product quality and customer service. Our highly trained staff will help you make the best choices for your needs while providing you with friendly service and expert advice.

Come see us today!

HOURS

MONDAY	8:30 - 6
TUESDAY	8:30 - 6
WEDNESDAY	8:30 - 6
THURSDAY	8:30 - 6
FRIDAY	8:30 - 6
SATURDAY	10 - 4
SUNDAY	CLOSED

Rental Equipment

Knee Walker  \$100/month	Wheelchair  \$100/month	Portable Oxygen  \$250/week \$500/month	Stationary Oxygen  \$200/month
CPAP  \$200/month	Lift Chair  \$300/month	Scooter  \$300/month	Hospital Bed  \$300/month

*ALL RENTAL EQUIPMENT REQUIRES A DEPOSIT

RENT WITH **YOU CAN HOME MEDICAL** 



Oxygen & Respiratory

Inogen One G5 System



G5 BACKPACK

G5 SINGLE BATTERY



One single battery · Concentrator · Standard 3-year warranty
AC/DC power supplies · Carry bag, user manual · Meets FAA requirements
Intelligent Delivery Technology® · Bluetooth System · Inogen Connect App

Comparable to OxyGo Next

Inogen One G4 System



One single battery · Concentrator · Standard 3-year warranty
AC/DC power supplies · Carry bag · Meets FAA requirements
Intelligent Delivery Technology® · Bluetooth System · Inogen Connect App

Comparable to OxyGo Fit

OxyGo Next



OxyGo NEXT is ground-breaking size and battery life have never been seen before at this high of an oxygen dose. Designed for 24/7 use Only 4.7 lbs.* with up to 13 hours of battery life** No carts or tanks required
Bluetooth-enabled units

Zen-O™



Zen-O™ portable oxygen concentrator can deliver up to 2 liters of oxygen in either pulse or continuous flow.

Zen-O lite™



Zen-O lite™ portable oxygen concentrator can deliver up to 1050ml of oxygen in pulse mode and has 5 flow settings.

OxyGo Fit



The 3 setting Intelligent Pulse dose OxyGo FIT is the smallest and lightest member of the OxyGo family! OxyGo FIT is about half the size of OxyGo, and is small enough and quiet enough to be worn on a patient's hip. Weighing in at 2.8 lbs.* with up to a 5 hour battery life** the OxyGo FIT Keeps Going—Everywhere.

Designed for 24/7 use
Only 2.8 lbs.* with up to 5 hours of battery life**

Devilbiss Stationary Unit




Made in the USA
Exclusive DeVilbiss OSD® (Oxygen Sensing Device) standard on every unit — increases accuracy and reliability for longer service intervals
Easy to assemble/disassemble two-piece cabinet for ease of maintenance
Accommodates up to a 50-foot cannula
Three Year Warranty

Inogen Stationary Unit



Weight: 18 lbs, Dimensions: 7 × 13 × 16.5 in
Oxygen Flow, Continuous Flow
Five flow settings from 1-5 liters per minute
Power Settings
AC Power Supply,
100-240V, 275W,
50-60Hz
Warranty: 3 Years or Lifetime



CPAP & Sleep Therapy

DreamStation Go

DreamStation Go is a ultra portable positive airway pressure (PAP) therapy system designed for the frequent traveler. Patients can rely on its unique blend of connectivity, pack-ability and powerful, convenient features.



AirMini™



- Smallest portable CPAP on the market today
- Effective waterless humidification for portable convenience
- Convenient control with the AirMini app
- Premium travel bag and accessories for therapy on the go

SoClean 2

You can use the SoClean 2 without needing to take equipment apart every day. Does not require water in order to ensure your equipment remains completely dry.

- New sleek design · Smaller footprint · Front display
- Indicator lights throughout the cycle
- Easy filter cartridge system · 2 Year Warranty
- SoClean 2 User Manual



AirSense™ 10 AutoSet™



- User-friendly controls and an intuitive interface make it simple to navigate settings and nightly sleep report
- AutoRamp™ feature delivers a low pressure to help you fall asleep, then ramps up to the required pressure

NUWAVE® CPAP Sanitizer

Automatically Kills 99.9% of Germs

The NUWAVE® CPAP Sanitizer is a revolutionary CPAP sanitizer which uses the proven platform of ozone gas to sanitize your CPAP equipment. NUWAVE® is different than all other CPAP sanitizers. NUWAVE® has been designed to be small as well as portable.



Staying Up To Date On The Newest Masks

RESMED Airfit F30



RESMED AirTOUCH F20



RESMED AIRFIT N30i



RESMED Airfit N20



Quattro™ Air



RESMED P10



AirMini™



Cleaning And Maintenance Accessories

Mask Wipes



Mask Cleaning Spray



Prescription required for CPAP Machines
NONE for Mask Purchase

Walkers & Wheelchairs



Nitro Aluminum Rollator



- Adjustable contoured backrest
- Cross-brace design for added safety
- Clip hold and side-by-side folding

Carbon Fiber Nitro



- Luxurious super lightweight carbon fiber frame
- Weighs only 12 lbs. for powerful performance with swift and effortless mobility
- Handle height easily adjusts with unique push button.

Nitro Zoom 22



Get the perfect fit and features with The Zoom. The Zoom has the best features with 8" wheels, large padded seat, and under seat pouch. Available in several colors.

Features: Locking Hand Brakes Large Padded Seat Under Seat Pouch Flip-Up Back Foldable



Nova Phoenix Rise Up



- Weight Capacity / Wt Cap: 300 lbs
- Weight 20.5 lbs · Overall Dimensions: 25" w x 33" d x 49.25"-53.25" h · Overall Dimensions When Folded: 12" w x 33" d x 49.25"-53.25" h
- Wheel Size: 10"(Front), 8"(Back)
- Approximate User Height: 5'2" – 6'0"

Stander Trust Walker



- It is light and elegant, both in use and when it is folded.
- It has a cross folding which makes it easy to fold and it takes little space and can stand for itself when it is folded

Lightweight Transport Chairs & Wheelchairs



- Comes with fixed desk length padded arms, quick release fold-down back, adjustable belt, durable thick nylon upholstery and back close-able pocket, removable and adjustable footrests, 8" front/rear wheels, rear wheel locks

Aloha Walker Bag	Black Walker Bag	Mobility Handbag	Mobility Handbag	Universal Mobility Bag	Mobility Bag	Mobility Bag
4001AP	4001BK	4006BB	4006BK	4009BK	4008GA	4008LP
Hanging Walker Pouch	Deluxe Cup Holder	Adjustable Cup Holder	Cane Holder	Cane Holder	Flashlight	Cane Hook
P4200SB	CH-1000	CH-1100	CH-4000R	CH-4100R	FL-2000	Cane Hook


Lift Chairs & Recliners



VivaLift!® - Atlas Plus



Color Options:



VivaLift!® - Tranquil



VivaLIFT!
POWER RECLINERS
by Pride Mobility Products Corp.



Power Headrest



Power Lumbar



USB Remote



Battery Backup



Footrest Extension



Power Backrest

Color Options:



VivaLift!® - Escape



Color Options:



VivaLift!® - Sierra



Color Options:



VivaLift!® - Elegance



Color Options:



Oasis Collection - LC-580iM



Color Options:




- Easy-to-use memory button hand control with USB charger port
- Premium furniture grade frames for added strength
- Sinuous no-sag back springs with the right amount of flexibility
- Removable filling allows for custom lumbar support
- Specially designed arms afford a luxurious feel and firm support

Heritage Collection - LC-358PW



Color Options:



- Industry-first hand control with USB charger port
- Standard head & arm covers
- Removable back
- Quiet & smooth lift system
- Wooden frame composed of engineered furniture grade laminate/hardwood for added strength and durability
- Wide range of sizes

Pride
Mobility Products Corp.

Lift Chairs & Recliners



Twilight- PR514-MLA



This unique fabric will enhance any relaxing experience with its one-of-a-kind characteristics:

- Breathable technology allows airflow to maintain constant temperature
- Plush and luxurious feel enhances comfort
- 300,000 Double Rubs – Developed to withstand above average wear and tear and constructed to be 10 times more durable than the industry standard for fabric
- Cleans effortlessly with soap and water, releases light stains with alcohol based household cleaner and disinfects with a 5:1 water/bleach solution for ease of care
- Innovative layering technology to repel moisture and liquids

Cloud



MaxiComforter



IMPERIAL



Easy Sleeper



Space Saver



Pub



DID YOU KNOW?

Although not considered a disability, people suffering from arthritis can benefit greatly from having a lift chair. Any slight amount of pressure on certain joints can lead to immense pain and stress. Lift chairs can reduce the amount of pressure and stress on the knees, back, and elbows, keeping arthritis from acting up when attempting to get up or sit down.

Lift chairs also greatly reduce swelling to your legs and feet.

ASK A REP ABOUT



- Total Control
- Position Control
- Twilight Control
- Custom Control



Powered Mobility Scooters & Chairs



Pride Go-Chair



Pride
Mobility Products Corp.

- 300 Lb. weight capacity
- Number one for portability
- Convenient under seat storage keeps your personal items close at hand
- Feather-touch disassembly easily disassembles into 5 lightweight pieces for transport and storage
- Your choice of six colors
- Heaviest piece is only 36 lbs.
- 18" x 17", 5-point swivel seat with latch release
- Speeds up to 3.7 mph



Wrangler

The Wrangler All-Terrain Scooter by Pride can take you anywhere! This heavy-duty scooter features two 3-horsepower motors, an unbeatable suspension system, and batteries that will last all day long. No other scooter on the market offers this much power for this price.



Pride
Mobility Products Corp.

- Travel Range: Up to 18.6 Miles Per Charge
- Seat: Ultra Padded Captains Chair
- Top Speed: 11.4 MPH
- Weight Capacity: 350 LBS
- Top Performing Outdoor Scooter
- Full Front / Rear LED Lighting Package

Go-Go LX



Pride
Mobility Products Corp.

Jazzy Passport



Pride
Mobility Products Corp.



Jazzy Zero Turn



Pride
Mobility Products Corp.

Grip Phone Holder



Weather Breaker



Rear-view Mirror



Safety Flag



Rear Basket



SCOOTER RENTALS

Great for vacations, outdoor, recovery, & so much more!

(303) 759-9525



Powered Mobility Scooters & Chairs



Folding Buzzaround




Air, sea, or train travel
300 LB weight capacity
Folds flat 12" / For Storage under Bed or furniture
Heaviest piece 47.5 LBS
In Stock : Blue & Red
Free Hard Shell Carry Bag
Requires minimal storage space
4- point carry handles for easy transport
12" height when folded (without seat)
USB port and water bottle holder on filler
15 AH lithium ion battery pack (standard) for a range of 18 miles
FREE CARRY ON- Hard Shell Rear Bag and Cup holder



iRide



This portable 3-wheel scooter features an easy-to-read digital odometer to keep track of your daily trips, manual handbrake with supplemental footbrake for added safety, 250-lb. weight capacity and speeds up to 3.7 mph and a twist-grip throttle. So, get out there and enjoy your life with this portable, lightweight, stylish iRide portable electric scooter!



LifeRider Envy



The World's lightest, least expensive, most maneuverable cash and carry powered wheelchair.

Unique, under seat storage basket standard
Tubular frame design
Colors: Cherry Red or Blue
Comfortable, full size stadium style seat
300 lb. weight capacity
Two 22Ah batteries (Now Standard)
Operating Range: Up to 15.5 miles
Ground Clearance: 2.5
Turning Radius: 26
Top speed: 3.5 mph
Designed exclusively for Golden, imported product



EZ Access Folding Ramps

Single-fold Ramp Hinge:

The Single-fold ramp features a continuous, no pinch, low profile hinge that runs the full length of the ramp to add strength and stability.

Ergonomically-designed Handles:


These make your SUITCASE Single-fold Ramp easy to carry, wherever you may go. Flexible and non-breakable, these handles offer ultimate comfort when carrying your portable ramp.

Slip-resistant Tread:

The extruded (built-in) tread allows mobility devices to retain traction in all weather conditions.

Self-adjusting Bottom Transition Plates:

Easily convert a mobility device from the SUITCASE Single-fold Ramp to any terrain with two independently self-adjusting bottom transition plates.



Harmar Vehicle Lifts

The AL100 Universal Scooter Lift is today's most popular outside scooter lift.

Features

- Sets the standard for all other Scooter Lifts
- Carries virtually all scooters
- Stars -N- Stripes Aluminum Deck
- Automatic hold down foot
- No scooter modifications
- Accommodates scooters with a wheelbase of up to 42"
- License plate mount included
- Swing-away option available





Seat Cushions & Pillows

FLIP Pillow	Wedge Cushion	Nova Bed Wedge
		
Elevating Leg Wedge	SitzRight Seat Cushion	Nova Cushion
		

Wellness & Pain

Trigger Point Massager	Blood Pressure Monitor	Finger Pulse Oximeters
		


You Can Home Medical offers hundreds of products to help with injury and pain.



CALL US FOR DETAILS

Aids To Daily Living

Handy Handle		Handy Bar		Great Grips	
					
Sugar Canes		Folding Canes		Standard Canes	
					
Universal Crutches	Nova Hip Kit	24" Reacher	32" Reacher	Pedal Exerciser	Hand Flex Star
					
Shoe Horn	Dressing Stick	Elastic Shoelaces	Sock Aid	Cast Protector	Heavy Duty Gait Belt
					





Medi Comfort Line



Compression Level	Length
15-20 mmHG	standard
20-30 mmHG	petite
30-40 mmHG	

REJUVA - In a world full of so many choices, Rejuva® provides high-quality, effective compression with stylized flair.



COMPRESSION ACCESSORIES

Butler

Simple step-in device to facilitate compression stocking application.



Why Compression?



varicose veins



diabetes



pregnant



airplane travel



swollen legs



surgery



stand for work



athletes



bed ridden



Universal Cervical Collar	Universal Shoulder Sling	CrissCross Sacro-Lumbar Support


Elbow Sleeve	Composite Wrist Brace	Knee Support

CEP Knee Sleeve	Hinged Knee Brace
Walker Boot	Stabilizing Ankle Brace

You Can offers orthopedic support products for your entire body.



We also carry a variety of hot/cold therapy products for reduced swelling and rapid healing.



Don't see the product you need? GIVE US A CALL!



Bemis Toilet Seat Riser



SEAT SECURITY: Patented Snap 2 Secure system is engineered to generate and maintain perfect clamping force to properly secure your seat so it will not wiggle or loosen and the self-sustaining hinges allows for easy cleaning, keeps seat up when in use and prevents ring slamming down






SHIELD PROTECTION: Innovative, patented "reverse funnel" design of the shield keeps everything in the bowl avoiding contact with waste

SAFETY & COMFORT: Seat height is increased/elevated by 3" to ease sitting and standing without bending for those with limited mobility from back, hip or knee issues

COMFORT & SECURITY: 21" width between the arms fits in virtually all bathrooms while providing extra width for the user and the heavy duty, solid plastic seat tested to 1,000 lbs and the ring's designed to eliminate painful pressure point contact

EASY TO USE: Arm width allows the seat ring to be raised in the upright position for easy cleaning and standing use of the toilet

Extra Wide for Comfort | Seat Holds 1,000 lbs

Bath/Shower Chair	Sliding Transfer Bench	Swivel Bath Seat	Bath Mat	Tub Grab Bar
Waterproof. Perfect for sitting in during showering or bathing.	Great for getting in and out of raised baths and showers with ease.	Swivels 360 degrees and locks in four 90 degree positions, which allows patients to be transferred into seat and swiveled into showering position.	Increase grip and eliminate slipping on bath/shower floors.	Extra boost of support for entering or exiting the tub/shower.
				

Hand Held Shower	Shower Chair Basket	Hinged Toilet Seat Riser	Toilet Safety Rails	Luxury Grab Bars	Adjustable Stool
					

Safety Tips: Ditch your throw rugs. These slippery floor coverings are the No. 1 cause of household falls. If you just can't part with them in the bathroom, secure loose rugs with double-faced tape, tacks or slip-resistant backings.

Take your time. The more you rush, the more likely you are to fall.



Home Safety

Sure Stand Security Pole with Handles



EZ Stand-N-Go



Bath Seat w/arms/back



Bath Seat Cushion



Padded Transfer w/back



Bath Seat Basket



Safe Home Bed Rail



Suction Cup Hook



Folding Commode



Commode Liners



Fall Facts: Falls are the leading cause of fatal and non-fatal injuries for older Americans. Falls threaten seniors' safety, independence and generate enormous economic and personal cost. Every 11 seconds, an older adult is treated in the emergency room for a fall; every 19 minutes, an older adult dies from a fall.

Attends Premier®



THE PREMIUM OVERNIGHT PROTECTION YOU DESERVE.

The premium protection of Attends Premier is focused on the unique needs of users and caregivers looking for a product that will perform no matter the situation, day or night.

Ultra-Plush Design

ConfidenceCuff™
Elastics

Odor Control

Dry-Lock™
Containment Core

SmartCool
Breath-ability

Premier Briefs



Premier Underwear



Premier Underpad



Catheters

These products + SO MUCH MORE!



How to Care for your Catheter

One-time use catheters and reusable catheters are available. For reusable catheters, be sure to clean both the catheter and the area where it enters the body with soap and water to reduce the risk of a UTI. One-time use catheters come in sterile packaging, so only your body needs cleaning before inserting the catheter.

You should also drink plenty of water to keep your urine clear or only slightly yellow. This will help prevent infection.

Empty the drainage bag used to collect the urine at least every eight hours and whenever the bag is full. Use a plastic squirt bottle containing a mixture of vinegar and water or bleach and water to clean the drainage bag.

Catheters



Attends Aloe Wipes



Reusable Underpad





Alterra 1385 Hi-Low Long-Term Hospital Bed

- Top-of-the-line parallel electronics system provides exceptionally smooth and quiet operation
- Synchronized high-low motors ensure the bed deck remains level even after frequent motion and long-term use
- Unique grid/rib deck prevents structure from sagging for durability
- Folding foot lock swings away when not in use to reduce tripping
- Shown with optional soft-touch assist rails and custom head and footboards for your convenience



Value Flex



Hi-Low



Premier



Basic Homecare Hospital Bed



- Headspring supports recessed to resist pinching and injuries; foot assembly recessed to reduce chance of pinching
- Light weight and 4 smooth-rolling casters (2 locking, 2 nonlocking) for easy transport
- DC motor box snaps off quickly for easy delivery, setup, and routine maintenance
- 9V battery backup for moving to the lowest position during power failures (batteries not included)
- Advanced remote with large buttons for easy patient use

Memory Foam




Innerspring/ Memory Foam



Gel Memory Foam



Latex



Low Profile



Deluxe Innerspring





YOU CAN HOME MEDICAL

4875 E Evans Ave STE 101, Denver, CO 80222 | 303-759-9525 | youcanhomemedical.com