

|                                               |                                                                                                           |
|-----------------------------------------------|-----------------------------------------------------------------------------------------------------------|
| Model                                         | Go Go® Sport <sup>1</sup>                                                                                 |
| Part Number                                   | 3-wheel: S73    4-wheel: S74                                                                              |
| Weight Capacity                               | 325 lbs.                                                                                                  |
| Maximum Speed <sup>2</sup>                    | Up to 4.7 mph                                                                                             |
| Ground Clearance <sup>4</sup>                 | 1.5" at motor / 0.75" at anti-tips                                                                        |
| Turning Radius <sup>4</sup>                   | 3-wheel: 36"    4-wheel: 41.5"                                                                            |
| Overall Length <sup>4</sup>                   | 3-wheel: 42"    4-wheel: 43"                                                                              |
| Overall Width <sup>4</sup>                    | 22.25"                                                                                                    |
| Seat-to-Ground Height                         | 20.5" - 22"                                                                                               |
| Front Tires                                   | 3-wheel: 3" x 8" solid<br>4-wheel: 3" x 9" solid                                                          |
| Rear Tires                                    | 3" x 9" solid                                                                                             |
| Range Per Charge <sup>1,2</sup>               | 3-wheel: Up to 12.75 miles<br>4-wheel: Up to 12 miles                                                     |
| Total Weight With Batteries                   | 3-wheel: 122.25 lbs.<br>4-wheel: 129.5 lbs.                                                               |
| Total Weight Without Batteries                | 3-wheel: 92.75 lbs.<br>4-wheel: 100 lbs.                                                                  |
| Heaviest Piece When Disassembled <sup>4</sup> | 3-wheel: 35.75 lbs. (front section)<br>4-wheel: 44 lbs. (front section)                                   |
| Standard Seating                              | Type: Compact, foldable<br>Weight: 23 lbs.<br>Material: Black vinyl<br>Dimensions: Width: 17", Depth: 17" |
| Motor Output                                  | 24-volt DC motor,<br>sealed transaxle                                                                     |
| Brakes                                        | Regenerative and<br>electromechanical                                                                     |
| Battery Requirements <sup>5</sup>             | (Two) 18 AH, 12 V, deep-cycle<br>Weight: 14 lbs. each                                                     |
| Battery Charger                               | Off-board, 2 amp                                                                                          |

## HCPCS Code: K0801 (Group 1, Heavy Duty)

- <sup>1</sup> Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- <sup>2</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- <sup>3</sup> Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- <sup>4</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- <sup>5</sup> Battery weight may vary based on manufacturer.

**PRIDE**  
MOBILITY

(US) 800-800-8586 • (Canada) 888-570-1113

pridemobility.com |    

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

041025

**GO GO**  
*SPORT*

**TRAVEL READY,  
NOW WITH TIGHTER  
TURNING!**

FEATURING EXCLUSIVE PATENTED **EZ TURN**  
TECHNOLOGY



**LIVE YOUR BEST®**

pridemobility.com |    

# NOW WITH TIGHTER TURNING!

FEATURING EXCLUSIVE PATENTED



Compact, foldable  
seat with armrests

Feather-touch  
disassembly



Delta tiller with  
charger port

High-intensity  
LED headlight

Interchangeable  
shroud panels  
in red and blue

The Go Go® Sport delivers high-performance and convenience on the go. With a 325 lb. weight capacity, a charger port in the tiller, standard front LED lighting and a maximum speed up to 4.7 mph., the Go Go Sport is feature rich and travel ready.

## Additional Features

- Pride's exclusive EZ Turn technology (4W unit only) provides 3-wheel turning ability while maintaining 4-wheel stability
- Convenient dual charging port to charge batteries on and off unit
- Deck basket for additional storage (4-wheel only)
- Front frame-mounted seat post offers maximum stability
- Front tiller-mounted basket
- Long-lasting, high-intensity, front LED headlight offers optimal pathway illumination
- Microprocessor-based controller offers optimal power management and added safety features
- Pride's exclusive black, non-scuffing tires

**Colors** (both color panel sets included)



Red

Blue

## Accessories\*

- |                               |                                                                          |
|-------------------------------|--------------------------------------------------------------------------|
| • Additional 18AH battery box | • Safety flag                                                            |
| • Cane/crutch holder          | • Under seat storage (assembly required)                                 |
| • Cell phone holder           | • Walker holder                                                          |
| • Cup holder                  | • Weather cover                                                          |
| • 50" Lap belt                | • Wishbone crutch holder                                                 |
| • Oxygen tank holder          | • XLR USB charger                                                        |
| • Rear basket                 | • Forearm crutch holder                                                  |
| • Rearview mirror(s)          | <b>Cup holder, rear basket, and saddlebag can be ordered as a bundle</b> |
| • Saddlebag                   |                                                                          |

## Military Patches



\* Only one rear accessory can be used at a time.

4-wheel model shown;  
3-wheel model available



Two sets of easily changeable,  
red and blue colored panels.

Feather-touch disassembly for  
quick transport or storage.

Exclusive EZ Turn Technology  
on 4-Wheel Model