Part Number	SC443E
Weight Capacity	300 lbs.
Maximum Speed[2]	Up to 4.23 mph
Ground Clearance[4]	1.29” at motor, 3.5” under deck
Turning Radius[4]	39.2”
Overall Length[4]	41.5”
Overall Width[4,6]	19.45”
Seat-to-Deck Height	16.38 – 17.83”
Seat-to-Ground Height	20.75” - 22.25”
Tires	Front: 2”x 7” Rear: 2.5”x 8”
Range[2,3]	12 Ah: Up to 7.1 miles at 300 lbs. Up to 8.9 miles at 220 lbs. 18 Ah: Up to 11.3 miles at 300 lbs. Up to 14.1 miles at 220 lbs.
Product Weight	83.6 lbs. without batteries 103 lbs. with 12 Ah batteries 111 lbs. with 18 Ah batteries
Heaviest Piece[4]	36.2 lbs. (front section)
Standard Seating	Type: Foldable, mesh back Material: Black vinyl Dimensions: 17”x17” Weight: 23 lbs.
Motor Output	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Brakes	Regenerative & electromechanical
Batteries	(Two) 12-volt, sealed lead-acid, non-spillable deep-cycle
Battery Weight[5]	12 Ah: 20 lbs. (standard) 18 Ah: 27.4 lbs.
Battery Charger	2 amp off-board charger
Warranty	Lifetime limited on frame; 2-year limited on Electronics; 2-year limited on Drivetrain Batteries: 6 months

HCPCS Code: K0800

1 Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
2 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
3 Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
4 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
5 Battery weight may vary based on manufacturer.
6 Dimensions listed are for scooter only and do not include seats or arms. Width may vary depending on seat options and accessories.

Featuring a new deck design with more leg room!



FDA Class II Medical Device[1]

GO-GO® ELITE TRAVELLER GEN 3



FDA Class II Medical Device[1]

Shown in Standard Cardinal Red

PRIDE® MOBILITY

pridemobility.com | 800.800.8586

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or – 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

pridemobility.com

082626

GO GO® ELITE TRAVELLER GEN 3

Featuring a new deck design with more leg room!



FDA Class II Medical Device[1]

FEATURING EXCLUSIVE PATENTED EZ TURN TECHNOLOGY



Convenient charging port in the tiller.



Intuitive control panel framed by a soft-touch delta tiller handle.



Tiller angle adjustment knob.

Go Go Elite Traveller® Gen 3 features EZ Turn technology and an ergonomically designed deck for improved leg room and pass-through space. A breathable mesh-back seat, wraparound delta tiller, large front basket and feather-touch disassembly also enhance comfort and convenience.

FEATURES

- **NEW** sloped deck design for improved comfort and **leg room**
- **Front and rear pathway lights**
- **EZ Turn** technology provides unrivaled tight turning without sacrificing the stability of a 4-wheel scooter
- Feather-touch **disassembly** allows for easy separation into 5 lightweight pieces for transport or storage
- Large front basket included
- Interchangeable shroud panels can be replaced with accessory colors shown below

INCLUDED COLOR PANELS (MATTE)

Cardinal Red

ACCESSORY COLOR PANELS (MATTE)

Black | Iceberg Blue | Robin's Egg Blue | Sugar Plum | Chell-o

ACCESSORIES*

- 2.1 amp XLR USB-A Charger
- 2.1 amp XLR USB-C Charger
- Cell phone holder
- Cup holder
- Lap belt
- Military patches
- Oxygen holder
- Rear basket
- Rear view mirror
- Safety flag
- Saddle bag
- Single/double rear mount 3"cane/crutch holder
- Under seat storage
- Walker holder
- Weather cover

* Only one rear accessory can be used at a time.